

# Laboratory diagnosis of pulmonary tuberculosis



**Dr. Kasumu Chisompola**

**MBS 240**

**12-07-2024**

# Objectives

- Know the causative agent of pulmonary TB
- Understand specimen collection and processing
- Understand the Laboratory methods for detection/diagnosis
- Understand recording and reporting of results

**What is pulmonary tuberculosis?**

- Tuberculosis is a bacterial infection caused by *Mycobacterium tuberculosis*
  - And other species of MTB complex including *M. bovis*, *M. africanum*, *M. canettii*
    - also known as tubercle bacilli or AFB.
  - Infection can be pulmonary or extrapulmonary
    - Can occur concurrently .
- Spread via airborne transmission from person to person in confined environment.

# Aerosol Formation: Spread of droplets from infected persons




- Coughing
- Singing
- Sneezing
- Talking

What sets *Mycobacterium tuberculosis* apart from other bacteria?

# The Mycobacterial Cell Wall

- Waxy cell wall composed of mycolic acid
- Outer coating serves as an efficient barrier and can endure unfavorable conditions
  - Disinfectants
  - Many drugs (antibiotics)
  - NaOH treatment during specimen processing
  - Staining procedures e.g. gram stain

# Mycobacterial cell wall





# Factors contributing to high tuberculosis infection rates

- High incidence of TB driven by:
  - HIV/AIDS
  - Poverty
  - Mis-use of anti-tuberculosis drugs
  - Development of drug resistant TB
  - Inadequate TB infection prevention and control (IPC) measures

# Impact of poverty on TB

- Crowded living conditions facilitate transmission
- Lack of education about symptoms
- Lack of infection control in homes and health care facilities
- Dependence on traditional healers
- Lack of access to health care sites
- Cost of health care
  - TB diagnosis & treatment is free in Zambia!

# Multidrug resistant tuberculosis (MDR-TB)

- First line drugs used for treatment of TB
  - Isoniazid, rifampicin, pyrazinamide & ethambutol
- MDR-TB is resistance to the most potent 1<sup>st</sup> line drugs
  - Isoniazid & rifampicin
- Treatment after failure with 2<sup>nd</sup> line drugs
  - Extensively drug resistant (XDR-)TB


# Specimen collection and processing

# Patients must be instructed properly

- Rinse mouth with clean water (boiled water).
- Open container but do not touch inside container or cap.
- Take 3–4 deep breaths, holding breath for 3-5 seconds before exhaling.
- Cough after the last exhalation.
- Empty sputum into container.
- Return specimens promptly to the lab.

# Containers for Sputum Collection

- Strong, unbreakable
- Leak proof, screw-capped with a water-tight seal
- Sterile
- Single-use
- Up to 50 ml capacity
- Translucent or clear material
- Easily-labeled walls



## Specimen containers must be labelled clearly

- Patients name, date of collection
- Complete request form to accompany specimen

## Specimens must be stored adequately

- Specimens must be promptly transported to lab
- Processed as soon as possible upon arrival at lab
- Stored in refrigerator (4°C) if not processed immediately

# Transporting specimens

- Specimens must be transported in a containment system
  - Primary containment
    - Collection container with screw cap top
  - Secondary containment
    - Specimen container in a sealable, biohazard bag
    - Place requisition in outside pouch of biohazard bag
  - Tertiary containment
    - Specimens in bags are placed in transport box (Styrofoam with fibreboard, plastic, or metal)
    - cold chain transport and keep specimens protected from light






# Quality Assessment of Sputum Specimens

- Characteristics of a good sputum specimen
  - Muroid or mucopurulent appearance
  - Minimum amounts of saliva
  - Optimal volume: 5ml-10ml
  - Minimum volume: 0.5 ml
- \*Record characteristics on form

# Assessing Quality of Sputum Specimens



# Personal protective equipment (PPE) MUST be worn



**N95 masks must be worn**



**Water resistant gowns and  
double gloves must be worn**

# Specimen processing

- All specimen processing should be carried out in a bio-safety cabinet (level 3)
  - Specimens should only be considered safe upon fixing onto a slide.
- Pre-treatment of sputum necessary before culture
  - NALC-NaOH commonly used method.

# Principles of processing

- Sputum specimens are viscous materials contaminated with normal flora
- Processing involves pre-treatment of the sputum specimens
  - Digestion: to free the TB bacilli from the mucus, cells or tissue in which they may be embedded
  - Decontamination: to eradicate normal flora that grow more rapidly than MTB and would interfere with the ability to recover MTB
  - Homogenization of the digested materials
  - Concentration of the TB bacilli by centrifugation before smear preparation and media inoculation


# Digestion-Decontamination

(N-Acetyl-L-Cysteine- Sodium Hydroxide method)


- NaOH
  - Decontaminating agent
- Na Citrate
  - Binds the heavy metal ions that might be present in the specimen that could inactivate the NALC
- N-Acetyl-L-Cysteine (NALC)
  - Digestion agent
  - Mucolytic agent aids in liquefying sputum to release trapped bacteria
  - more efficient decontamination of normal flora and releases AFB

# AFB-smear microscopy

- Performed on sputum
- Can be direct or concentrated sputum
- Most mycobacteria are Acid fast bacilli (AFB)
  - Resistant to decolourisation during staining
- Microscopy usually using light microscope
  - Ziehl-Neelsen stain
- Primary means of diagnosis in resource limited settings
- Although ZN microscopy has been phased out and replaced with FM



**ZN stain of AFB:**  
AFB retain red dye after  
Decolourising step



**AFB under light microscope**  
Following ZN staining




# Grading of ZN smear

Number of AFB	Fields	Report
No AFB	In 100 immersion fields	No AFB seen
1-9	In 100 immersion fields	Exact No
10-99	In 100 immersion fields	1+
1-10	Per-field examination	2+
>10	Per-field examination	3+

# Fluorescence microscopy

- Requires smears to be stained with fluorescent stain
  - Auramine-rhodamin
- More sensitive than ZN microscopy
  - Will detect small numbers of bacilli
- More slides can be examined compared to ZN
- Lower resolution required than ZN

# Fluorescence LED microscopy of *Mycobacterium tuberculosis*



# LED fluorescence Microscopy

## TB grading

- No AFB: report as “No AFB seen”
- 1-19 AFB/40 fields: actual number recorded
- 20-199 AFB/40 fields: 1+
- 5-50 per field: 2+
- >50 per field: 3+

# Culture

- Gold standard for diagnosis
  - more sensitive
- Antigen detection strips specific for MTBC used following culture e.g. Capilia

## 1. Solid culture

- Can take up to 6-8 weeks for growth to be observed

## 2. Liquid culture-BACTEC MGIT system

- Automated culture system
- Can take up to 20 days for growth
- Expensive and only available at referral level

# Solid culture

- Lowenstein-Jensen (LJ)  
for isolation of MTB
- Contains glycerol,  
potato flour, salts &  
eggs
  - Glycerol replaced with  
Na pyruvate for *M. bovis*
- Malachite green inhibits  
gram +ve bacteria



# Liquid culture

- Performed in automated culture system
  - Mycobacteria Growth Indicator Tube (MGIT)
- BACTEC MGIT system provides rapid detection for TB growth
- Mycobacteria use  $O_2$  during growth
  - $O_2$  released from indicator in tube bottom
  - Fluorescence indicates growth of bacteria
- Incubates samples at  $37^\circ\text{C}$

# Length of Incubation in MGIT

- Protocol length for growth detection can be from 1 to 56 days
  - Default protocol length is 42 days for detection of TB from processed sputum specimens
  - Positive specimens ~20 days
- Specimens with no growth after 42 days can be reported as negative
  - Subcultures or AFB stains are not required unless flakes of growth are seen in tube





MGIT tubes




BACTEC MGIT machine

- Advantages of culture
  - Detects small numbers of organisms
    - Can detect as little as 10 bacilli/ml
    - Confirms diagnosis of TB in HIV patients
  - Allows species identification
  - Allows for drug sensitivity testing
- Limitations of culture
  - Slow growth, hence long turn-around time
  - Expensive
  - Limited laboratories
    - 3 reference laboratories providing liquid cultures

# Antigen detection test for culture

- AFB smear positive actively growing culture
  - Liquid
  - Solid
- Interpretation after 15 minutes



# Interpretation


## Negative

No *Mycobacterium tuberculosis* antigen



## Positive

*Mycobacterium tuberculosis* complex positive



## Invalid

It is recommended that the specimen should be re-tested



## Antigen Detection Test

# Drug susceptibility testing (DST)

- Important for testing sensitivity of TB to anti-TB drugs
- Colonies from solid and liquid culture can be used for testing
- Sensitivity testing carried out at referral level using automated BACTEC MGIT system
- Strains resistant to the first line drugs isoniazid and rifampicin are considered to be MDR

# Molecular methods

- Provide definitive identification of MTBC
- Amplifies and detects genomic material
- Rapid detection
- More sensitive than microscopy and culture
- Limitations:
  - Expensive equipment and reagents
  - Requires trained personnel to operate
- Systems for identification include:
  - Line probe assay: identification of different MTBC
  - GeneXpert: RIF resistant marker indicating MDR
  - PCR

# GeneXpert

- Automated diagnostic test for ID of MTB
- Cartridge based PCR system
- Detects rifampicin resistance
- Highly sensitive & specific
  - Recommended by WHO for use in HIV-TB co-infection





# Reporting of results


- Laboratory results must be returned to the clinician promptly
- Results should indicate:
  - If MTBC isolated
  - Any other microorganisms isolated
  - Drug susceptibility test results
- Specimen rejection must be reported with reasons
  - Details on container not matching request form
  - missing request form
  - no patient details on specimen

# Why is the TB laboratory important?

- Patient and community care
  - Allows for the patient to be placed on appropriate treatment
  - Shortens infectious period
- Diagnosis of TB
  - Allows for timely and accurate diagnosis
  - Patient management for drug therapy
  - Monitoring and treatment of drug resistant TB

In summary

# Isolation and Identification of Mycobacteria



# Further reading

- **Murray, Medical Microbiology** 6<sup>th</sup> edition.
- **Mims' Medical Microbiology** 4<sup>th</sup> Ed